

Alkaline Phosphatase (ALP)

Interpretive Summary

Description: Alkaline Phosphatase (ALP) is primarily an indicator of cholestatic liver disease. It also increases with severe bone destruction and due to steroid induction.

Decreased ALP

Common Causes

- Not clinically significant
- Artifact
 - Hemolysis
 - Assay performed on EDTA plasma

Increased ALP

Common Causes

- Cholestasis, either intrahepatic or extrahepatic (dogs >cats or horses)
 - Impaired bile flow due to:
 - Hepatic necrosis
 - Hepatocellular swelling
 - Cushing's disease (dogs)
 - Hepatic lipidosis (cats)
 - Diabetes mellitus
 - Gall bladder mucocoele/stones
 - Hyperthyroidism (cats)
 - Neoplasia
 - Inflammatory
 - Cholangiohepatitis/cholecystitis
 - Pancreatitis
- Induction by drugs or hormones
 - Corticosteroids
 - Phenobarbital
 - Thyroxine
- Increased osteoblastic activity in bone
 - Young animals (increases up to 3-4 times level expected in mature dogs)
 - Osteosarcoma, other bone neoplasia

Uncommon Causes

- Drugs (dogs)
 - Phenytoin
 - Primidone
- Toxic
 - Ragwort (horses)
 - Alsike clover (horses)
 - Mycotoxin
- Increased osteoblastic activity in bone (dogs)
 - Fracture repair
 - Rickets
 - Hyperparathyroidism

- Benign familial hyperphosphatasemia (Siberian Huskies)
- Benign hyperphosphatasemia (Scottish Terriers)
 - Adults >6 yo, may have higher ALP values than other breeds when matched for age
- Miscellaneous
 - Neoplasia (benign and mammary tumors in dogs)
 - Post-suckling pups and foals
 - Increased placental ALP in late term pregnancy (cats)

Related Findings

- Cholestasis
 - Increased ALT, AST, GGT, bilirubin
 - Abdominal radiographs may show:
 - Enlarged liver
 - Gall stones or calcification of gall bladder (rare)
 - Abdominal ultrasound may show:
 - Enlarged, hyperechoic liver
 - Abnormal gall bladder
 - Enlarged/obstructed
 - Mucocoele
 - Stones in the gall bladder
 - Histopathology/cytology findings consistent with cholestatic liver disease
- Cushing's disease
 - Decreased urine specific gravity
 - Stress leukogram: increased neutrophils and monocytes, decreased lymphocytes and/or eosinophils
 - Adrenal function tests consistent with Cushing's disease
- Hepatic lipidosis
 - ALP:GGT ratio increased compared to other liver diseases in cats
- Diabetes mellitus
 - Increased serum glucose and glucosuria
 - Increased fructosamine
- Hyperthyroidism
 - Increased T4, free T4, free T4 by equilibrium dialysis
- Neoplasia
 - Enlarged/irregular liver on radiographs and/or ultrasound
 - Cytology/histopathology findings consistent with neoplasia
- Pancreatitis
 - Increased amylase and lipase
 - Increased Spec cPL® or Spec fPL®
- Osteosarcoma
 - Proliferative bone lesions on radiographs
 - Cytology/histopathology findings consistent with neoplasia

Additional Information

Physiology

- Alkaline phosphatase refers to a group of enzymes that catalyze hydrolysis of phosphate esters in an alkaline environment in vitro.
- Cellular function is poorly described, but these are membrane-bound enzymes present in most tissues.
- High activities are present in liver, bone, intestine, kidney, and placenta.
- Changes in the serum level of ALP are attributed primarily to hepatic (L-ALP) and bone (B-ALP) isoenzymes since intestinal, renal, and placental isoenzymes have extremely short half lives.

- Hepatic ALP is found mainly in liver canalicular cell membranes and increases with biliary disease, especially with cholestasis.
- In dogs, glucocorticoids (endogenous or exogenous) induce an isoenzyme in the liver (C-ALP) that contributes to serum increases
- In cats, the hepatic isoenzyme has such a short half life that even minor increases are considered significant.
- In horses, serum ALP levels are highly variable, and are of limited diagnostic use

References

- Latimer KS, Mahaffey EA, Prasse KW, eds. *Duncan and Prasse's Veterinary Laboratory Medicine: Clinical Pathology*, 4th ed. Ames, IA: Blackwell; 2003.
- Gallagher AE, Panciera DL, Panciera RJ. Hyperphosphatasemia in Scottish terriers: 7 cases. *JVIM* 2006; 20(2): 418-421.
- Nestor DN, Holan KM, Johnson CA et al. Serum alkaline phosphatase activity in Scottish Terriers versus dogs of other breeds. *JAVMA* 2006; 228(2): 222-228
- Stockham SL, Scott MA. *Fundamentals of Veterinary Clinical Pathology*, 2nd ed. Ames, IA: Blackwell; 2008.
- Willard MD, Tvedten H, eds. *Small Animal Clinical Diagnosis by Laboratory Methods*, 4th ed. St. Louis, MO: Saunders; 2004.

Last Updated 11/1/2013